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Referenze

1. Andersson M, Alvarez-Cermeno J, Bernardi G, Cogato I, Fredman P, Frederiksen J et al (1994) Cerebrospinal fluid in the diagnosis of multiple sclerosis: a consensus report. *J Neurol Neurosurg Psychiatry* 57:897–902
2. Sellebjerg F, Christiansen M (1996) Qualitative assessment of intrathecal IgG synthesis by isoelectric focusing and immunodetection: interlaboratory reproducibility and interobserver agreement. *Scand J Clin Lab Invest* 56:135–143
3. Nespolo A, Bernardi G (1995) Lo studio delle proteine liquorali: aspetti fisiopatologici e nuove prospettive diagnostiche. Corso CEFAR 1995: Le proteine dal laboratorio alla clinica
4. Reiber H (1995) External quality assessment in clinical neurochemistry: survey of analysis for cerebrospinal fluid (CSF) proteins based on CSF/serum quotients. *Clin Chem* 41:256–263
5. Teunissen CE, Malekzadeh A, Leurs C, Bridel C, Killestein J (2015) Body fluid biomarkers for multiple sclerosis—the long road to clinical application. *Nat Rev Neurol* 11:585–596
6. Granerod J, Ambrose HE, Davies NW, Clewley JP, Walsh AL, Morgan D et al (2010) Causes of encephalitis and differences in their clinical presentations in England: a multicentre, populationbased prospective study. *Lancet Infect Dis* 10:835–844
7. Gable MS, Sheriff H, Dalmau J, Tilley DH, Glaser CA (2012) The frequency of autoimmune N-methyl-D-aspartate receptor encephalitis surpasses that of individual viral etiologies in young individuals enrolled in the California Encephalitis Project. *Clin Infect Dis* 54: 899–904
8. Gresa-Arribas N, Titulaer MJ, Torrents A, Aguilar E, McCracken L, Leypoldt F et al (2014) Antibody titres at diagnosis and during follow-up of anti-NMDA receptor encephalitis: a retrospective study. *Lancet Neurol* 13:167–177
9. Freedman MS, Thompson EJ, Deisenhammer F, Giovannoni G, Grimsley G, Keir G et al (2005) Recommended standard of cerebrospinal fluid analysis in the diagnosis of multiple sclerosis: a consensus statement. *Arch Neurol* 62:865–870
10. Polman CH, Reingold SC, Banwell B, Clanet M, Cohen JA, Filippi M et al (2011) Diagnostic criteria for multiple sclerosis: 2010 revisions to the McDonald criteria. *Ann Neurol* 69:292–302
11. Reiber H, Felgenhauer K (1987) Protein transfer at the blood cerebrospinal fluid barrier and the quantitation of the humoral immune response within the central nervous system. *Clin Chim Acta* 163: 319–328
12. Hutchinson M (2013) CSF oligoclonal bands are important in the diagnosis of multiple sclerosis, unreasonably downplayed by the McDonald criteria 2010: commentary. *Mult Scler* 19:719–720
13. Masjuan J, Alvarez-Cermeno JC, Garcia-Barragan N, DiazSanchez M, Espino M, Sadaba MC et al (2006) Clinically isolated syndromes: a new oligoclonal band test accurately predicts conversion to MS. *Neurology* 66:576–578
14. Zipoli V, Hakiki B, Portaccio E, Lolli F, Siracusa G, Giannini M et al (2009) The contribution of cerebrospinal fluid oligoclonal bands to the early diagnosis of multiple sclerosis. *Mult Scler* 15: 472–478
15. Milo R, Miller A (2014) Revised diagnostic criteria of multiple sclerosis. *Autoimmun Rev* 13:518–524
16. Thompson AJ, Banwell BL, Barkhof F, Carroll WM, Coetzee T, Comi G, et al. (2017) Diagnosis of multiple sclerosis: 2017 revisions of the McDonald criteria. *Lancet Neurol* 17:162-173.
17. Joseph FG, Hirst CL, Pickersgill TP, Ben-Shlomo Y, Robertson NP, Scolding NJ (2009) CSF oligoclonal band status informs prognosis in multiple sclerosis: a case control study of 100 patients. *J Neurol Neurosurg Psychiatry* 80:292–296
18. Siritho S, Freedman MS (2009) The prognostic significance of cerebrospinal fluid in multiple sclerosis. *J Neurol Sci* 279:21–25

19. Thangarajh M, Gomez-Rial J, Hedstrom AK, Hillert J, AlvarezCermeno JC, Masterman T et al (2008) Lipid-specific immunoglobulin M in CSF predicts adverse long-term outcome in multiple sclerosis. *Mult Scler* 14:1208–1213
20. Garcia-Barragan N, Villar LM, Espino M, Sadaba MC, GonzalezPorque P, Alvarez-Cermeno JC (2009) Multiple sclerosis patients with anti-lipid oligoclonal IgM show early favourable response to immunomodulatory treatment. *Eur J Neurol* 16:380–385
21. Villar LM, Costa-Frossard L, Masterman T, Fernandez O, Montalban X, Casanova B et al (2015) Lipid-specific immunoglobulin M bands in cerebrospinal fluid are associated with a reduced risk of developing progressive multifocal leukoencephalopathy during treatment with natalizumab. *Ann Neurol* 77:447–457
22. Mandrioli J, Sola P, Bedin R, Gambini M, Merelli E (2008) A multifactorial prognostic index in multiple sclerosis. Cerebrospinal fluid IgM oligoclonal bands and clinical features to predict the evolution of the disease. *J Neurol* 255:1023–1031
23. Gastaldi M, Zardini E, Franciotta D (2017) An update on the use of cerebrospinal fluid analysis as a diagnostic tool in multiple sclerosis. *Expert Rev Mol Diagn* 17:31–46
24. Majed M, Fryer JP, McKeon A, Lennon VA, Pittock SJ (2016) Clinical utility of testing AQP4-IgG in CSF: guidance for physicians. *Neurol Neuroimmunol Neuroinflamm* 3:e231
25. Wingerchuk D, Banwell B, Bennett J, Cabre P, Carroll W, Chitnis T et al (2015) International consensus diagnostic criteria for neuromyelitis optica spectrum disorders. *Neurology* 85:177–189
26. Stoeck K, Sanchez-Juan P, Gawinecka J, Green A, Ladogana A, Pocchiari M et al (2012) Cerebrospinal fluid biomarker supported diagnosis of Creutzfeldt-Jakob disease and rapid dementias: a longitudinal multicentre study over 10 years. *Brain* 135:3051–3061
27. Fokke C, van den Berg B, Drenthen J, Walgaard C, van Doorn PA, Jacobs BC (2014) Diagnosis of Guillain-Barré syndrome and validation of Brighton criteria. *Brain* 137:33–43
28. Joint Task Force of the EFNS and the PNS (2010) European Federation of Neurological Societies/Peripheral Nerve Society guideline on management of chronic inflammatory demyelinating polyradiculoneuropathy: report of a joint task force of the European Federation of Neurological Societies and the Peripheral Nerve Society—first revision. *J Peripher Nerv Syst* 15:1–9
- 29 32. Martino G, Grimaldi LM, Muiola L, Filippi M, Martinelli V, Comi G et al (1990) Discontinuous distribution of IgG oligoclonal bands in cerebrospinal fluid from multiple sclerosis patients. *J Neuroimmunol* 30:129–134
30. Bertolotto A, Malentacchi M, Capobianco M, di Sapio A, Malucchi S, Motuzova Y et al (2016) The use of the 25 Spotte needle markedly reduces post-dural puncture headache in routine neurological practice. *Cephalalgia* 36:131–138
31. Cruickshank AM (2001) ACP best practice no 166: CSF spectrophotometry in the diagnosis of subarachnoid haemorrhage. *J Clin Pathol* 54:827–830
32. Steele RW, Marmer DJ, O'Brien MD, Tyson ST, Steele CR (1986) Leukocyte survival in cerebrospinal fluid. *J Clin Microbiol* 23: 965–966
33. Tibbling G, Link H, Ohman S (1977) Principles of albumin and IgG analyses in neurological disorders. I. Establishment of reference values. *Scand J Clin Lab Invest* 37:385–390
34. Thompson E (1988) *The CSF proteins: a biochemical approach*. Elsevier, Amsterdam
35. Gregoire SM, van Pesch V, Goffette S, Peeters A, Sindic CJ (2006) Polymerase chain reaction analysis and oligoclonal antibody in the cerebrospinal fluid from 34 patients with varicella-zoster virus infection of the nervous system. *J Neurol Neurosurg Psychiatry* 77: 938–942

36. Koedel U, Fingerle V, Pfister HW (2015) Lyme neuroborreliosis epidemiology, diagnosis and management. *Nat Rev Neurol* 11:446–456
37. Davis LE, DeBiasi R, Goade DE, Haaland KY, Harrington JA, Harnar JB et al (2006) West Nile virus neuroinvasive disease. *Ann Neurol* 60:286–300
38. Hottenrott T, Dersch R, Berger B, Rauer S, Eckenweiler M, Huzly D et al (2015) The intrathecal, polyspecific antiviral immune response in neurosarcoidosis, acute disseminated encephalomyelitis and autoimmune encephalitis compared to multiple sclerosis in a tertiary hospital cohort. *Fluids Barriers CNS* 12:27
39. Reiber H, Lange P (1991) Quantification of virus-specific antibodies in cerebrospinal fluid and serum: sensitive and specific detection of antibody synthesis in brain. *Clin Chem* 37:1153–1160
40. Dorries R, Ter Meulen V (1984) Detection and identification of virus-specific, oligoclonal IgG in unconcentrated cerebrospinal fluid by immunoblot technique. *J Neuroimmunol* 7:77–89
41. Olsson T, Kostulas V, Link H (1984) Improved detection of oligoclonal IgG in cerebrospinal fluid by isoelectric focusing in agarose, double-antibody peroxidase labeling, and avidin-biotin amplification. *Clin Chem* 30:1246–1249
42. Davies G, Keir G, Thompson EJ, Giovannoni G (2003) The clinical significance of an intrathecal monoclonal immunoglobulin band: a follow-up study. *Neurology* 60:1163–1166
43. Franciotta D, Bergamaschi R, Amato MP, Zardini E, Persico A, Portaccio E et al (2005) Clinical correlations of CSF single IgG bands. *J Neurol* 252:1274–1275
44. Ferraro D, Franciotta D, Bedin R, Solaro C, Cocco E, Santangelo M et al (2017) A multicenter study on the diagnostic significance of a single cerebrospinal fluid IgG band. *J Neurol* 264:973–978
45. Franciotta D, Avolio C, Lolli F (2005) Between-laboratory variability in oligoclonal IgG band numbering. *Clin Chem* 51:270–272
46. Hegen H, Auer M, Zeileis A, Deisenhammer F (2016) Upper reference limits for cerebrospinal fluid total protein and albumin quotient based on a large cohort of control patients: implications for increased clinical specificity. *Clin Chem Lab Med* 54:285–292