## **RCR 2020**

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#### **DEMENZE**

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### III. BIBLIOGRAFIA

- Betsou F, Gaignaux A, Ammerlaan W, Norris PJ, Stone M. Biospecimen Science of Blood for Peripheral Blood Mononuclear Cell (PBMC) Functional Applications. *Current Pathobiology Reports* 2019;7:17–27. doi: 10.1007/s40139-019-00192-8.
- Böyum A. Separation of leukocytes from blood and bone marrow. Introduction. *Scand J Clin Lab Invest Suppl* 1968;97:7.
- Chen H, Schürch CM, Noble K, et al. Functional comparison of PBMCs isolated by Cell Preparation Tubes (CPT) vs. Lymphoprep Tubes. *BMC Immunol* 2020;21:15. doi: 10.1186/s12865-020-00345-0.
- Corkum CP, Ings DP, Burgess C, Karwowska S, Kroll W, Michalak TI. Immune cell subsets and their gene expression profiles from human PBMC isolated by Vacutainer Cell Preparation Tube (CPT™) and standard density gradient. *BMC Immunol* 2015;16:48. doi: 10.1186/s12865-015-0113-0.
- Grievink HW, Luisman T, Kluft C, Moerland M, Malone KE. Comparison of Three Isolation Techniques for Human Peripheral Blood Mononuclear Cells: Cell Recovery and Viability, Population Composition, and Cell Functionality. *Biopreserv Biobank* 2016;14:410-415. doi: 10.1089/bio.2015.0104.
- Jiang CY, Niu Z, Green MD, et al. It's not 'just a tube of blood': principles of protocol development, sample collection, staffing and budget considerations for blood-based biomarkers in immunotherapy studies. *J Immunother Cancer* 2021;9:e003212. doi: 10.1136/jitc-2021-003212.
- Mallone R, Mannering SI, Brooks-Worrell BM, et al. Isolation and preservation of peripheral blood mononuclear cells for analysis of islet antigen-reactive T cell responses: position statement of the T-Cell Workshop Committee of the Immunology of Diabetes Society. *Exp Immunol* 2011;163:33-49. doi: 10.1111/j.1365-2249.2010.04272.x.
- Turner RJ, Geraghty NJ, Williams JG, et al. Comparison of peripheral blood mononuclear cell isolation techniques and the impact of cryopreservation on human lymphocytes expressing CD39 and CD73. *Purinergic Signal* 2020;16:389-401. doi: 10.1007/s11302-020-09714-1.
- Ulmer AJ, Scholz W, Ernst M, Brandt E, Flad HD. Isolation and Subfractionation of Human Peripheral Blood Mononuclear Cells (PBMC) by Density Gradient Centrifugation on Percoll. *Immunobiology* 1984;166:238-250. doi: 10.1016/S0171-2985(84)80042-X.
- van Beem RT, Hirsch A, Lommerse IM, et al. Recovery and functional activity of mononuclear bone marrow and peripheral blood cells after different cell isolation protocols used in clinical trials for cell therapy after acute myocardial infarction. *EuroIntervention* 2008;4:133–138. doi: 10.4244/eijv4i1a21.